PRODUCT INFORMATION SHEET

Intended Use

TransFix is a solution that preserves cellular antigen expression (cluster of differentiation markers), and cell morphology within biological samples such as: peripheral whole blood and cerebrospinal fluid. Samples treated with TransFix can be maintained for up to 14 days prior to analysis by flow cytometry.

TransFix is a General Diagnostic Reagent.

Summary and Principles

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes in a blood sample. Any delay in testing the sample, such as transport from the collection site to analysis location, can have a negative impact on results. Sample preservation provides a solution which addresses delays between sample collection and testing.

Subsets of leukocytes can be distinguished by cell surface antigens using fluorescent antibodies and flow cytometry. Qualitative and quantitative changes in leukocyte subsets are used to identify and monitor immunodeficiency and hematological diseases [1]. The TransFix stabilizer acts by preserving the cell surface antigens of lymphocyte subsets until processing and analysis can be performed.

TransFix is available in bulk format: supplied in 1ml and 20ml aliquots in polypropylene vials, or as a 'sample storage tube': consisting of a 1.2ml polypropylene vial containing 0.2ml TransFix.

Reagents

TransFix is a clear green liquid containing formaldehyde and other chemicals.

Disposal

TransFix contains formaldehyde which should be disposed of in accordance with local regulations. Avoid disposing into drainage systems and the environment.

Once the TransFix 'sample storage tube' contains a biological sample it must be disposed in accordance with local regulations regarding clinical waste.

Precautions and Warnings

1. TransFix is intended for use as specified in this document and is a General Diagnostic Reagent, IVD Class I.
2. For stabilizing peripheral blood samples, anti-coagulated blood must be kept at room temperature prior to treatment with TransFix.
3. Incubation times or temperatures other than those specified will lead to erroneous results.
4. Treat the sample with TransFix as quickly as possible to ensure the best results and no later than 6 hours after collection.
5. Do not use TransFix after the expiration date on the tubes and packaging.
6. Do not dilute or add other components to TransFix.
7. Do not use cell viability stains on samples treated with TransFix as they are fixed instantaneously.
8. Do not re-use TransFix Sample Storage Tubes.
9. TransFix treated samples and all materials coming into contact with it should be handled as if capable of transmitting infection.
10. Avoid contact of TransFix treated samples with the skin and mucus membranes. The cell preservative is an irritant and any contact should be washed off with soap and water immediately.
11. TransFix does not contain any antimicrobial reagents. Microbial contamination should be avoided, or erroneous results may occur.
12. SDS can be obtained at www.cytomark.co.uk or by calling +44(0)1280 827460.

Indications of Product Deterioration

1. Cloudiness or precipitate visible in the TransFix.
2. Color change of TransFix from a clear green solution.
3. TransFix change from liquid to solid.
4. If indications of product deterioration occur, do not use the TransFix and contact Cytomark immediately on: +44(0)1280 827460 or support@cytomark.co.uk.

Storage Conditions and Stability

All TransFix products are shipped at ambient temperature (2 - 25°C). Additional insulation may be required for shipment during extreme temperature conditions. TransFix products are supplied in sealed foil pouches.

Unused **TransFix Sample Storage Tubes** are stable at 2 - 25°C for up to 12 months or until the expiration date on the label. *N.B.* 2 - 8°C storage is advisable where possible.

Unused **TransFix bulk** is stable at 2 - 25°C for up to 24 months or until the expiration date on the label. *N.B.* 2 - 8°C storage is advisable where possible.

Instructions for Use

1. Collect blood by venipuncture into an EDTA vacuum tube according to normal phlebotomy practice.
2. For **TransFix Sample Storage Tubes**, mix the blood sample, remove 1ml and add to the tube.
3. For **Bulk TransFix**, carefully remove the blood collection tube cap and measure the volume of anti-coagulated whole blood within the vacuum tube. Pipette into

the blood collection tube the appropriate volume of TransFix at the ratio of 0.2ml TransFix per 1ml of blood.

4. Replace the cap on the blood collection tube, ensuring that there is no leakage.
5. Mix the treated sample carefully. Do not vortex.
6. Store / transport the TransFix treated blood for up to 14 days at 2 - 8°C, **or**, Store / transport the TransFix treated blood at 18-25°C for up to 3 days, and at 2-8°C for a up to a further 4 days.
7. Incubate the treated sample at room temperature (18-25°C) for 15 minutes before use.
8. Re-suspend the cells thoroughly by **gently inverting the tube at least 10 times**. Heavier cells and blood components will sediment over time, forming two distinct layers; this is normal.
9. **Take care when opening blood collection tubes as there may be a build-up of blood in the caps.**
10. Perform analysis by flow cytometry in accordance with the manufacturer's instructions. A 'stain, lyse-no wash' sample preparation method is recommended. Analyze the treated sample within 6 hours before being returned to 2 - 8°C storage for future use, if necessary.
11. **When treated with Transfix the dilution factor must be accounted for when calculating absolute cell counts. Adjust the absolute cell count by multiplying the output by 1.2.**

Notes:

1. Use caution when implementing automatic gating strategies as light scatter positions of cells stabilized by TransFix may differ from those of untreated cells.
2. Studies have shown that increased levels of hemolysis, icterus and lipemia do not affect the results.

TransFix is quality tested using a panel of markers including CD3, CD4, CD8 and CD45. A certificate of analysis and a certificate of conformity can be provided with every batch of TransFix and is available at www.cytomark.co.uk.

It is recommended that all antibody conjugates are validated in association with TransFix prior to use. A list of antibodies compatible with TransFix can be found at www.cytomark.co.uk.

References

1. Evaluation of stabilized blood cell products as candidate preparations for quality assessment programs for CD4 T-cell counting. Bergeron et al, Clinical Cytometry, Vol. 50, 2002, 86-91.

Ordering Information

Please call Cytomark on +44(0)1280 827460 or email support@cytomark.co.uk for assistance. Additional information can be found online at www.cytomark.co.uk.

| Product Description | Catalogue Numbers |
|--|-------------------|
| TransFix Bulk (1 x 1ml tube) | TFB-01-1 |
| TransFix Bulk (10 x 1ml tube) | TFB-01-10 |
| TransFix Bulk (50 x 1ml tube) | TFB-01-50 |
| TransFix Bulk (1 x 20ml tube) | TFB-20-1 |
| TransFix Sample Storage Tube (2 x 1.2ml tube) | TF-01-2 |
| TransFix Sample Storage Tube (10 x 1.2ml tube) | TF-01-10 |
| TransFix Sample Storage Tube (25 x 1.2ml tube) | TF-01-25 |
| TransFix Sample Storage Tube (50 x 1.2ml tube) | TF-01-50 |

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Cytomark, une division de Caltag Medsystems Ltd.
Caltag Medsystems Ltd, Buckingham, MK18 1TF, UK